




GALIX

BIOMEDICAL INSTRUMENTATION, INC.

ECG-PS Computerized ECG and Stress Test System

- Very High Resolution 12-lead ECG Data Acquisition module
- Specialized, Powerful Measuring Tool
- Unique menu driven report edition Tool
- More than 150 parameters are automatically measured on each beat



Computerized ECG

- Full disclosure module for all ECG data stored of any lead
- Ventricular Late Potential Analyzer Software option
- Export/Import to/from other systems
- Unlimited database



Computerized Stress Test

- Control capabilities of different treadmills and bicycles
- Pre-programmed standard protocols + user friendly tool for easy protocols creation
- Study can be reviewed at any time while performing test using review screen; we do not miss any single beat
- Unique Report Generator; customized report can be issued within seconds.
- Full Print Preview adds flexibility and reliability
- Several drop-down menus with pre-programmed phrases speeds up report edition
- Wandering baseline, muscle and line noises digital filters allow very clean
- ECG waveforms even on highest stages of Stress Study

2555 Collins Ave., Suite C-5, Miami Beach, FL 33140, USA. Phone: 305-534-5905. Fax: 305-534-8222

sales@galix-gbi.com

www.galix-gbi.com